

Inhalation Sciences





DissolvIt® – evaluation report and presentation of data
from the FDA funded study on Inhalation Sciences *in vitro* module

Maria Malmlöf, Per Gerde, Manoush Masarrat

2025-05-28

Introduction

- The DissolvIt® system is an *in vitro* test model/dissolution apparatus that is built to resemble relevant lung physiology for evaluation of the dissolution- and absorption of inhalable drugs.
- Drug particles are evenly distributed over glass cover slips and dissolved into a mucus simulant and then absorbed into flowing perfusate, thereby creating *in vivo* like conditions and generating time-concentration curves with T_{\max} and C_{\max} as measurable parameters.



Project aims

- The aims of the FDA funded project were to:
 - **Aim 1:** Evaluate the discriminatory ability of the DissolvIt system using different formulations with known **differences or similarities**.
 - **Aim 2:** Directly compare DissolvIt data to IPL data in rat *ex vivo* as well as to clinical data *in vivo*.
 - **Aim 3:** Investigate the potential for *in vivo* predictability of DissolvIt data by performing physiologically-based biopharmaceutical modeling (PBBM) .

Inhalation products tested in Dissolv^{It} within this project

16 test products

6 manufactured micronized
dry powders (pure API)



9 commercially
available DPIs




1 commercial available
pMDI



Methods: Aerosolization and dose deposition (PreciseInhale®) dissolution testing (DissolvIt®)



6 manufactured micronized
dry powders (pure API)



9 commercially available
DPIs






1 commercial available
pMDI

Test product

5



Deposition/exposure flow
Actuation flow



Dissolution testing in DissolvIt


**Aerosolization and dose deposition onto glass cover
slips with PreciseInhale**

Aim 1: Evaluate the discriminatory ability of the DissolvIt system using different formulations with known differences or similarities

ISAB code	Test product	Evaluations performed
IS194	Budesonide, DISCUS, 3.29 μm (X50)	Data evaluation of APSD and dissolution of products where the API PSD and API manufacturing method are varied.
IS195	Budesonide, DISCUS, 2.14 μm (X50)	
IS196	Budesonide, UMAX, 1.54 μm (X50)	

ISAB code	Test product	MMAD (μm)	GSD
IS194	Budesonide, DISCUS, 3.29 μm (X50)	2.34 ± 0.00	2.33 ± 0.11
IS195	Budesonide, DISCUS, 2.14 μm (X50)	1.78 ± 0.18	2.14 ± 0.05

Aim 1: Budesonide – different sizes, differences expected to be seen




The smaller particles (IS195) are dissolved faster, that is demonstrated by identifying a shorter T_{max} , higher normalized C_{max}

Dissolv*It* detects differences

*Statistically significant difference $p < 0.05$, Student's t-test, two-sided, assuming similar variance

Aim 1: Budesonide – different sizes, differences expected to be seen



....and higher values of the cumulative dissolution.

DissolvIt detects differences

*Statistically significant difference $p < 0.05$, Student's t-test, two-sided, assuming similar variance

ISAB code	Test product	Normalized C_{max} (%/mL)	T_{max} (min)	Cumulative dissolution at 15 min
IS194	Budesonide, DISCUS, 3.29 μm (X50)	8.3 ± 0.2	3.0 ± 0.0	40.4 ± 0.7
IS195	Budesonide, DISCUS, 2.14 μm (X50)	$9.8^* \pm 0.4$	$2.0^* \pm 0.0$	$48.1^* \pm 2.5$

Aim 1: Evaluate the discriminatory ability of the DissolvIt system using different formulations with known differences or similarities


ISAB code	Test product	Evaluations performed
IS194	Budesonide, DISCUS, 3.29 μm (X50)	Data evaluation of APSD and dissolution of products where the API PSD and API manufacturing method are varied.
IS195	Budesonide, DISCUS, 2.14 μm (X50)	
IS196	Budesonide, UMAX, 1.54 μm (X50)	

Aim 1: Budesonide – different manufacturing methods

ISAB code	Test product	Manufacturing method	MMAD (μm)	GSD
IS195	Budesonide, 2.14 μm (X50)	DISCUS	1.78 ± 0.18	2.14 ± 0.05
IS196	Budesonide, 1.54 μm (X50)	UMAX	1.59 ± 0.12	1.98 ± 0.13

- DISCUS = Dispersive crystallization with ultrasound, gives more "normal" crystalline particles
- UMAX = ultrasound mediated amorphous to crystalline transition, gives crystalline particles with round morphology, more rugosity

Aim 1: Budesonide – different manufacturing methods, differences expected to be seen



*Statistically significant difference $p < 0.05$, Student's t-test, two-sided, assuming similar variance

ISAB code	Normalized C_{max} (%/mL)	Cumulative dissolution and absorption (%) after								
		20 s	40 s	1 min	2 min	3 min	5 min	8 min	15 min	30 min
IS196, Bud UMAX 1.54 μm (X50)	14.6 \pm 1.9*	0.3 \pm 0.1	1.2 \pm 0.2*	2.7 \pm 0.4*	8.2 \pm 1.2*	14.3 \pm 2.0*	26.0 \pm 3.5*	40.4 \pm 5.3*	60.2 \pm 6.1*	76.4 \pm 4.5
IS195, Bud DISCUS 2.14 μm (X50)	9.8 \pm 0.4	0.2 \pm 0.0	0.9 \pm 0.0	1.9 \pm 0.0	5.7 \pm 0.2	9.8 \pm 0.3	17.6 \pm 0.7	28.3 \pm 1.2	48.1 \pm 2.5	70.6 \pm 5.3

DissolvIt detects differences

The more rough UMAX particles were expected to dissolve faster.

Aim 1: Evaluate the discriminatory ability of the DissolvIt system using different formulations with known differences or similarities

ISAB code	Test product	Evaluations performed
IS406	Symbicort Turbohaler 320/9 (Bud/For F), DPI	APSD determination and dissolution testing in DissolvIt for both APIs in a brand name product versus a generic product.
IS407	Bufomix Easyhaler (320/9) (Bud/For F), DPI	
IS408	Seretide Evohaler FP/SX (250/25), pMDI	APSD determination of both APIs in an <i>in vitro</i> and an <i>ex vivo</i> set-up. Dissolution testing of both APIs in DissolvIt and generation of lung absorption data in IPL (<i>ex vivo</i>) for both APIs. Comparison of the generated <i>in vitro</i> and <i>ex vivo</i> data with existing <i>in vivo</i> data. PBB modeling of the DissolvIt generated FP data.



Test product	MMAD of Bud (μM) (n=3)	GSD of Bud (n=3)	MMAD of For (μM) (n=3)	GSD of For (n=3)
Symbicort Turbohaler, IS406	3.11 ± 0.17	1.77 ± 0.01	3.12 ± 0.19	1.76 ± 0.01
Bufomix Easyhaler, IS407	3.31 ± 0.15	1.77 ± 0.00	3.79 ± 0.22	1.63 ± 0.00



Brand name (Symbicort Turbohaler) and generic product (Bufomix Easyhaler) – similarities expected to be seen[#]



Symbicort Turbohaler is a brand name product in Europe, and Bufomix Easyhaler is now an approved generic product in Europe. Both contains two APIs, budesonide and formoterol.



No statistically significant differences (Student's t-test, two-sided assuming similar variance).

Dissolv/It confirms similarities seen in clinical data[#].


Aim 2: Directly compare DissolvIt data to IPL data in rat *ex vivo* as well as to clinical data *in vivo*

PreciseInhale

DissolvIt – *in vitro*



- Non-biological system physiologically resembling the lung
- Comprises an artificial air-blood barrier with a mucus simulant
- System perfused with a blood simulant; dissolved and absorbed API detected in the single-pass perfusate over time (4 h)



Clinical exposures – *in vivo*



- Regional lung dosing with PreciseInhale
- Blood samples collected for 24 h
- API analysis in plasma samples
- Clinical study performed: Gerde, P., et al., *Regional lung targeting with a fluticasone/salmeterol aerosol using a bolus breath hold method of the PreciseInhale® system: A first evaluation in humans.* Eur J Pharm Sci, 2024. **196**: p. 106742.

Isolated and Perfused Lung (IPL) – *ex vivo*






- Specific exposure of the rat lung
- Lung ventilated and perfused during the experiment (2 h)
- Perfusate analysis provides lung specific PK-data
- Substance remaining in lung completes a mass balance

Aim 2: Directly compare DissolvIt data to IPL data in rat *ex vivo* as well as to clinical data *in vivo*

ISAB code	Test product	Evaluations performed
IS406	Symbicort Turbohaler 320/9 (Bud/For F), DPI	APSD determination and dissolution testing in DissolvIt for both APIs in a brand name product versus a generic product.
IS407	Bufomix Easyhaler (320/9) (Bud/For F), DPI	
IS408	Seretide Evohaler FP/SX (250/25), pMDI	APSD determination of both APIs in an <i>in vitro</i> and an <i>ex vivo</i> set-up. Dissolution testing of both APIs in DissolvIt and generation of lung absorption data in IPL (<i>ex vivo</i>) for both APIs. Comparison of the generated <i>in vitro</i> and <i>ex vivo</i> data with existing <i>in vivo</i> data. PBB modeling of the DissolvIt generated FP data.

System	Test Product	IS #	MMAD FP (μm) GSD FP	MMAD Sal (μm) GSD Sal
DissolvIt – <i>in vitro</i>	Seretide Evohaler (250/25), pMDI	408	3.88 ± 0.10	3.97 ± 0.10
			1.71 ± 0.04	1.69 ± 0.05
IPL – <i>ex vivo</i>	Seretide Evohaler (250/25), pMDI	408	3.98 ± 0.20	4.07 ± 0.20
			1.91 ± 0.08	1.89 ± 0.08
Clinical – <i>in vivo</i>	Seretide Evohaler (250/25), pMDI	n/a*	4.22 ± 0.11	4.55 ± 0.02
			1.98 ± 0.02	2.22 ± 0.1



DissolvIt and IPL comparison – *in vitro ex vivo* correlation (IVEVC)



- Similar dissolution/absorption curves for both FP and Sal as well as same order of magnitude for C_{max} for both FP and Sal from Seretide Evohaler 250/25 in DissolvIt and IPL
- Similar dissolved and absorbed % API from Seretide Evohaler 250/25 in DissolvIt and IPL after 2 h for both FP (11% in DissolvIt, 14% in IPL) and Sal (75% in DissolvIt, 71% in IPL)
- The small differences seen can be explained by the system set-up/design such as time for dose deposition and lipid distribution in the air-blood barrier

Direct data comparison is possible
between DissolvIt and IPL

IVIVC – *in vitro/ex vivo in vivo* correlation



- Circulatory clearance of Sal and FP from Seretide Evohaler 250/25 pMDI in DissolvIt (green), IPL (orange) and clinic* (blue).
- DissolvIt data compares more directly to IPL data and describes dissolution and absorption of substance in the lung only whereas clinical data shows the additional effect of the substance passing into the blood circulation.
- Together with PBB modeling the DissolvIt FP data was used to predict a clinical profile.


Aim 3: Investigate the potential for *in vivo* predictability of DissolvIt data by performing physiologically-based biopharmaceutical modelling



- The FP dissolution data generated in DissolvIt for Seretide Evohaler 250/25 pMDI together with dose deposition data from literature were used for PBB modeling of a FP plasma concentration-time profile.
- This calculated profile was compared with the observed clinical data (Gerde, P., et al., *Regional lung targeting with a fluticasone/salmeterol aerosol using a bolus breath hold method of the PreciseInhale® system: A first evaluation in humans*. Eur J Pharm Sci, 2024. **196**: p. 106742.).

Can Dissolv/*t* predict clinical data?

PBB modelling of FP in Seretide Evohaler 250/25 pMDI



PBB = physiologically-based biopharmaceutical, FP = fluticasone propionate, Cp = plasma concentration,
Red dotted line (calculated) = Dissolv/*t* dissolution data after PBB modeling, Black dots (observed) = Observed clinical data

Dissolv/*t* dissolution method predicted a clinically relevant dissolution rate for FP in Seretide Evohaler 250/25 pMDI

Conclusions

- Dissolv/*t* detects expected differences in dissolution/absorption originating from
 - Different particle sizes (budesonide)
 - Difference in API manufacturing method (UMAX and DISCUS manufacturing of budesonide)
- Dissolv/*t* detects similarities in dissolution/absorption as expected for
 - Brand name and generic product (budesonide and formoterol from Symbicort Turbohaler and Bufomix Easyhaler)
- Dissolv/*t* generates concentration curves as well as cumulative absorption of FP and Sal in Seretide Evohaler (250/25) very similar to those generated in IPL.
- **Ample evidence that Dissolv/*t* can correctly detect potential differences in dissolution/absorption profiles originating in alterations of the test formulations. Dissolv/*t* also has the potential to generate data that can be used to predict clinical plasma profiles.**

Acknowledgement

- FDA
- pharm-analyt (LC-MS/MS analysis)
- Karolinska Institutet (SEM images)
- Emmace (development of Mimetikos Preludium module and PBB modeling)



Q/A